

## TECHNICAL INFORMATION FOR **ECO RECEIPT 04**

|                                 |  |
|---------------------------------|--|
| TYPE:                           | Acidic solution of surfactants   |
| DELIVERY<br>SHAPE:              | Single component water system  |
| FEATURES<br>AND USE:            | <p><b>ECO RECEIPT 04</b> is used to remove corrosion from the surface of metals and metal alloys. It is primarily used for cleaning rust from iron and steel.</p> <p><b>ECO RECEIPT 04</b> also presents one type of metal pre-protection. After drying of the means any standard color for metal can be applied on the treated surface.</p> <p>When is necessary to keep the original color of metal <b>RECEIPT ECO 04</b> is used along with <b>ECO IMPREGNIR M</b>, means for the protection of metals and metal alloys from the corrosion.</p> |
| PREPARATION FOR<br>APPLICATION: | Before the use <b>ECO RECEIPT 04</b> should be diluted with cold water in a volume ratio of 1:5, 1:4 or 1:3 (medium : water), depending on the amount of corrosion that needs to be removed.   |
| APPLICATION<br>METHOD:          | <b>ECO RECEIPT 04</b> is applied to the metal by spraying with pneumatic sprayer or coating with brush or sponge. After applying the means surfaces that are much corroded should be rinsed with water and wiped, then treated again.  |
| SAFEGUARD<br>MEASURES:          | <b>ECO RECEIPT 04</b> is acidic detergent, so it is necessary to use the appropriate safety equipment. In case the concentrated detergent comes in contact with exposed skin wash the spot with cold water. If it accidentally gets into the eyes flush them immediately with cold water. Rinsing with sodium bicarbonate solution is also recommended. Surfaces that are accidentally sprayed with this means should be washed with water. For professional use only!   |
| PACKAGING:                      | 10L Plastic cans and 1L plastic bottles.   |
| STORAGE:                        | <b>ECO RECEIPT 04</b> is stored outdoors and indoors at temperatures above 4°C.  |
| PRODUCT<br>LIFETIME:            | <b>ECO RECEIPT 04</b> can be kept in the original packaging for 12 (twelve) months.  |